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AR

Safety Information and Warnings for Users of the PAL System

Edition 7 – 2026 – Original Instructions

Contents:

This safety guide raises awareness of potential safety issues and general points for consideration for CTC Analytics representatives during installation and repair of the PAL System or parts of it (following the life cycle principle), as well as for the end user in the lab during the learning phase and in routine work.

It is important to read this section first before operating the PAL System.

Product Warranty and Disclaimer	2
General Considerations	3
Environmental Conditions	3
Electrical Hazards	4 - 5
Laser Safety Information	6
Other Hazards	6 - 8
Working with Toxic or Other Harmful Compounds	8
Biological Hazards	9
Maintenance	9 - 10
Disposal	10
Regulatory Compliance	11
FCC Compliance Statement	11
Safety Alerts and Important Information	12
Commonly Used Symbols for Instrument Markings	13 - 14

Product Warranty and Disclaimer

CTC Analytics AG reserves the right to make improvements and/or changes to the product(s) described at any time without prior notice.

CTC Analytics AG makes no warranty of any kind pertaining to this product, including but not limited to implied warranties of merchantability and suitability for a particular purpose.

Under no circumstances shall **CTC Analytics AG** be held liable for any coincidental damage or damages arising as a consequence of or from the use of this document.

© 2026 **CTC Analytics AG**. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine-readable form without prior written permission from **CTC Analytics AG**, except as permitted under copyright laws.

CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

General Considerations

The PAL System User Manual and related documents must be consulted by the user under all circumstances before a unit is put to use.

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The user should be aware that if the equipment is used in a manner not specified by the manufacturer, the protective and safety features of the equipment may be impaired.

The repair of instrument failures caused by operation in a manner not specified by the manufacturer is expressly excluded from the standard warranty and service contract coverage.

When, for technical reasons, it is necessary to work on instrument parts which may involve a potential hazard (moving parts, components under voltage, etc.), the authorized representative of CTC Analytics AG must be contacted. In general, this type of situation arises when access to the parts is only possible using a tool. Should an operator perform such a maintenance operation, he/she must have received proper training to carry out that specific task.

Environmental Conditions

The environmental conditions in which the PAL System can be operated are listed in the corresponding section of the PAL User Manual; see the section 'Specifications'.

Electrical Hazards

Every analytical instrument has specific hazards, so be sure to read and comply with the following precautions. They will help ensure the safe, long-term use of your PAL System.

The installation category (over-voltage category) for this instrument is Level II. The Level II category pertains to equipment receiving its electrical power from a local level, such as an electrical wall outlet.

Connect the PAL System only to instruments complying with IEC 61010 safety regulations. The power line and the connections between the PAL System and other instruments used in the configuration setup of the total analytical system must maintain good electrical grounding. Poor grounding represents a danger for the operator and may seriously affect instrument performance.

Do not connect the PAL System to power lines that supply devices of a heavy duty nature, such as motors, refrigerators and other devices that can generate electrical disturbances.

Use only fuses of the type and current rating specified. Do not use repaired fuses and do not short-circuit the fuse holder.

The supplied power cord must be inserted into a power outlet with a protective earth contact (ground). When using an extension cord, make sure that the cord also has an earth contact.

If the supplied power cord does not fit the local electrical socket and a replacement or adapter has to be purchased locally, make sure that only a certified power cord is used. Any power cord used should be certified by the appropriate local authorities.

Pay attention not to leave any cable connecting the PAL System and the chromatographic system or the power cord close to heated zones, such as the injector or detector heating blocks or the GC hot air vents.

Always replace any cable showing signs of damage with another one provided by the manufacturer. Safety regulations must be respected.

Do not change the external or internal grounding connections. Tampering with or disconnecting these connections could endanger you and/or damage the PAL System.

The instrument is properly grounded in accordance with these regulations when shipped. You do not need to make any changes to the electrical connections or the instrument's chassis to ensure safe operation.

The combination of a PAL System with an LC/MS System also requires that the safety measures described by the LC/MS System manufacturer be followed. Detailed instructions for the safety grounding on the LC/MS system are outlined in the corresponding operating/installation manual.

CTC Analytics recommends using a grounding cable connected on one side to the Injection Valve, Loop or any other suitable direct metallic contact, and on the other side to an appropriate ground point on the LC/MS System. This supplementary grounding measure will complement the safety strategy of the LC/MS System manufacturer.

Electrical Hazards

Do not turn the instrument on if you suspect that it has incurred any kind of electrical damage. Instead, disconnect the power cord and contact a CTC Analytics representative for a product evaluation. Do not attempt to use the instrument until it has been evaluated. Electrical damage may have occurred if the PAL System shows visible signs of damage, exposure to any liquids or has been transported under severe stress.

Damage can also result if the instrument is stored for prolonged periods under unfavorable conditions (e.g. subjected to heat, moisture, etc.).

Ensure that the power supply/controller unit is always placed in a clean and dry position. Avoid any liquid spills in the vicinity.

Before attempting any type of maintenance work, always disconnect the power cord(s) from the power supply(ies) if optional devices are installed. Capacitors inside the instrument may still be charged even though the instrument is turned off.

To avoid damaging electrical parts, do not disconnect an electrical assembly while power is being supplied to the PAL System. Once the power is turned off, wait approximately 30 seconds before you disconnect an assembly.

The instrument includes a number of integrated circuits. These circuits may be damaged if exposed to excessive line voltage fluctuations and/or power surges or electrostatic charges.

Never try to repair or replace any components of the instrument that are not described in the PAL User Manual without the assistance of a CTC Analytics representative.

There are no operator-serviceable or replaceable parts inside the power supply(ies) or in the PAL System.

If a power supply is not functioning, contact a CTC Analytics representative.

The power supplies for the PAL System and the Peltier Stack Module have the symbols I/O on the power switch label to indicate ON/OFF.

If a Peltier Stack Module is installed in combination with a PAL System, then a second power supply is active in the complete system. Turning OFF the two power supplies or pulling the two power cords in an emergency will stop the entire PAL System.

It is important that the power supply(ies) be in a location where the power ON | OFF switch is accessible and easy to operate, and where it is possible to unplug the AC power cord from the power supply/wall outlet in case of emergency.

Laser Safety Information

Safety Warning for Class 1 Laser Product

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Warning:

The installed laser device is a Class 1 Laser Product. Class 1 laser devices are not considered to be hazardous when used for their intended purpose. The following statement is required to comply with US and international regulations:

Caution:

Use of controls, adjustments or performance of procedures other than those specified herein may result in hazardous laser light exposures.

The selected Class 1 laser for the PAL module Barcode Reader complies with the following regulations:

- 21 CFR1040.10 and 1040.11 except for conformance with IEC 60825-1 Ed. 3., as described in Laser Notice No. 56, dated May 26, 2019
- EN60825-1:2014+ A11:2021
- IEC60825-1:2014/AMD1:2017

The software contains a built-in safety time limit such that the laser scanning mechanism cannot be operated in AIM mode for more than 5 continuous seconds

Other Hazards

To avoid injury and possible infection through contamination during PAL System operation, keep your hands away from the syringe.

Do not operate the PAL System without the safety guard. The safety guard must be installed for safe operation.

Do not place any objects inside the area of the safety guard. Keep away from the area around the safety guard during operation of the PAL System.

Other Hazards

Danger of crushing to fingers and hands. To avoid injury, keep your hands away from moving parts during operation. Turn off the power to the PAL System if you need to reach inside a mechanically powered system with moving parts.

Do not operate on the instrument components that form part of the work area of the PAL System when this is in motion.

To avoid injury, observe safe laboratory practice when handling solvents, changing tubing or operating the PAL System. Know the physical and chemical properties of the solvents you use. See the MSDS (Material Safety Data Sheets) from the manufacturer of the solvents being used.

When using the PAL System, follow the generally accepted procedures for quality control and method development.

When using the PAL System in the field of chromatographic analysis, if a change in the retention of a particular compound, in the resolution between two compounds or in peak shape is observed, immediately determine the reason(s) for the changes. Do not rely on the separation results until you determine the cause of a change.

Use caution when working with any polymer tubing under pressure:

- Always wear eye protection when near pressurized polymer tubing.
- Do not use polymer tubing that has been severely stressed or kinked.
- Do not use polymer tubing, in particular PEEK or Tefzel tubing, when using tetrahydrofuran (THF), dimethylsulfoxide (DMSO), chlorinated organic solvents, concentrated mineral acids such as nitric, phosphoric or sulfuric acids, or any related compounds.

Do not use vials without a sealing cap, or microtiter or deepwell plates without a plate seal. Vapor phase from organic solvents can be hazardous and flammable. Acidic vapor phase can cause corrosion to critical mechanical parts.

Good glass quality is essential whenever samples are heated and agitated in glass vials. High pressure may build up in closed vials depending on the type of solvent, application conditions (e.g. temperature settings) and possible chemical reactions. When heating liquids in vials, it is strongly recommended not to set the temperature higher than 10 °C above the boiling point of the solvent used. Note that even solid materials may contain volatile compounds and residual solvents such as water (in the form of moisture), which can build up excess vapor pressure.

Do not re-use headspace vials. When washing the vials, micro-cracks can form which will weaken the glass and increase the risk of vial breaking.

Other Hazards

In case of a single fault situation where the temperature control of the Agitator fails, there is the potential danger that the device will heat up in an uncontrolled manner until it reaches the cut-off temperature of the overtemperature fuse, in this case, 240 °C.

Based on this single fault scenario, when working with flammable solvents, the user must ensure that the solvent used has a fire point which is 25 °C higher than that of any of the heated surfaces that it could potentially come in contact with.

When filling-up a standard reservoir or replacing a solvent such as a washing solvent, remove the solvent reservoir bottle from the system to avoid a possible spill over the instrument. Depending on the physical, chemical or hazardous properties of the solvent, use the appropriate protective measures for handling.

Working with Toxic or Other Harmful Compounds

Before handling dangerous substances (toxic, harmful, etc.), please read the hazard indications and information reported in the Material Safety Data Sheet (MSDS) supplied by the manufacturer referring to the relevant CAS (Chemical Abstract Service) number. The PAL System requires the use of several chemical products with different hazard characteristics, which are present in vials and syringes. Before using these substances or replacing the syringe, please read the hazard indications and information reported in the MSDS supplied by the manufacturer referring to the relevant CAS number.

During sample preparation, please refer to local regulations for the ventilation conditions of the work room.

All waste materials must be collected and eliminated in compliance with the local regulations and directives in the country where the instrument is being used.

Biological Hazards

In laboratories where samples with potential biological hazards are handled, the user must label any equipment or parts thereof which may become contaminated with biohazardous material. The appropriate warning labels are included with the shipment of the instrument. It is the user's responsibility to label the relevant parts of the instrument.

When working with biohazardous materials, it is the user's responsibility to fulfill the following mandatory requirements:

- Instructions on how to safely handle biohazardous material must be provided.
- Operators must be trained and made aware of the potential dangers.
- Personal protective equipment must be provided.
- Instructions must be provided on what to do in case operators are exposed to aerosols or vapors during normal operation (within the intended use of the equipment) or in case of single fault situations such as a broken vial. The protective measures should consider potential contact with the skin, mouth, nose (respiratory organs) and eyes.
- Instructions for decontamination and safe disposal of the relevant parts must be provided.

It is the user's (operator's) responsibility to handle hazardous chemicals or biological compounds, including (but not limited to) bacterial or viral samples and the associated waste, safely and in accordance with international and local regulations.

Maintenance

Any external cleaning or maintenance must be performed with the PAL System turned off and the power cord disconnected. Avoid using solvents and spraying on electrical parts. For the removal of potentially dangerous substances (toxic, harmful, etc.), read the hazard indications and information reported in the MSDS (Material Safety Data Sheet) supplied by the manufacturer referring to the relevant CAS (Chemical Abstract Service) number. Use proper protective gloves.

When working with hazardous materials such as radioactive, biologically hazardous material, etc., it is important to train all operators how to respond in case of spills or contamination.

Depending on the class of hazardous material, the appropriate measures have to be taken immediately. Therefore, the chemicals or solvents needed for decontamination have to be on hand.

Maintenance

Any parts of the equipment which can potentially be contaminated, such as the sample vial rack, syringe tool, wash module, etc., should be cleaned regularly. The waste solvent from cleaning and any hardware which needs to be disposed of has to be properly eliminated with all the necessary precautions, abiding by national and international regulations.

When preparing for decontamination, ensure that the solvent or chemical to be used will not damage or react with the surface, dye (color) of the instrument, table or other nearby objects.

If in doubt, please contact your CTC Analytics representative to verify the compatibility of the type or composition of solvents with the PAL System.

It is the user's (operator's) responsibility to handle hazardous chemicals or biological compounds, including (but not limited to) bacterial or viral samples and the associated waste, safely and in accordance with international and local regulations.

Disposal

Do not dispose of this equipment or parts thereof unsorted in municipal waste. Follow local municipal waste regulations for proper disposal provisions to reduce the environmental impact of waste electrical and electronic equipment (WEEE).

European Union customers: Call your local customer service representative responsible for the PAL System for complimentary equipment pick-up and recycling.

Warning:

The customer has to ensure that the PAL System has not been contaminated by any hazardous chemical or biological compounds, including (but not limited to) bacteria or viruses.

Any part which has had direct contact with the analytical sample must be identified and must undergo an appropriate decontamination procedure prior to shipping for disposal. Potentially dangerous components are: Syringes, Valves, Loops, Tubing, Vials and Well Plates.

Any critical parts sent for disposal must be handled according to national laws for hazardous compounds.

The customer and the service engineer/technician are fully responsible for enforcing these requirements.

CTC Analytics AG will hold the representative and/or customer responsible if these regulations are not observed.

Regulatory Compliance

CTC Analytics performs complete testing and evaluation of its products to ensure full compliance with applicable domestic and international regulations.

CTC Analytics declares, under sole responsibility, that the product as originally delivered complies with the requirements of the following applicable European directives and carries the CE marking accordingly:

EMC Directive: 2014/30/EU
Machinery Directive: 2006/42/EC, Machinery Regulation: EU (2023/1230
RoHS Directive: 2011/65/EU, amended by Delegate Directive 2015/863/EU

and conforms with the following product standards:

Electromagnetic Compatibility (EMC):

IEC 61326-1:2020, EN 61326-1:2021
EN 55011:2016+A1:2017, +A11:2020+A2:2021, CISPR 11:2015 + A1:2016, CISPR11:2024
FCC part 15 Subpart B, §15.107(a); §15.109(a)

Safety:

IEC 61010-1:2010/AMD1:2016, EN 61010-1:2010 + A1:2019/AC:2019
UL 61010-1:2012 Ed.3.1., CSA C22.2 No. 61010-1:2012/A1:18

FCC Compliance Statement

CTC Analytics AG hereby declares that our product has been tested and complies with the requirements of FCC Part 15 Subpart B (Class B).

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Safety Alerts

Safety Alerts and Important Information

Make sure that you follow the precautionary notices presented in this guide and in any additional documentation such as the PAL User Manual. Safety and other special notices include the following:

Notice on the Proper Use of CTC Analytics Instruments.

In compliance with international regulations: Use of this instrument in a manner not specified by CTC Analytics AG could impair the safety and protective features of the instrument.

Warning:

For your safety, and in compliance with international regulations, the physical handling of this CTC Analytics instrument requires a team effort to lift and/or move the equipment. This instrument is too heavy and bulky for one person to handle safely. Before lifting or moving the PAL System, first secure all moving parts to avoid crushing hands or tipping over the unit.

Warning:

This is the general safety warning symbol and safety alert word to prevent actions that could cause personal injury. It highlights a general hazard or signifies the presence of another specific safety symbol. This is usually followed by an appropriate safety message. When you see a safety alert on your instrument or in the documentation, please carefully read the safety instructions before proceeding.

Note:

Note that highlights information necessary to prevent damage to software, loss of data, invalid test results or instrument damage; might contain information that is critical for optimal performance of the system or other important information about a task.

Commonly Used Symbols

Symbol	Description
	Caution, or Refer to User Manual
	Caution, Risk of Needle-stick Puncture
	Caution, Hot Surface or High Temperature
	Danger of Crushing to Fingers and Hands
	
	Laser Warning, Barcode Reader
	Biological Hazard
	Direct Current
	Alternating Current

Commonly Used Symbols

Symbol	Description
	Protective Conductor Terminal, Ground
	Fuse
	Electrical Power ON Used with Main PAL Power Supply
	Electrical Power OFF Used with Main PAL Power Supply
	Caution, Risk of Electrical Shock (High Voltage)
	Disposal: Do not dispose of in municipal waste. Follow local waste regulations to reduce electrical and electronic waste (WEEE).

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© 2026 **CTC Analytics AG**. Tutti i diritti riservati. La presente pubblicazione non può, né interamente né in parte, essere copiata, fotocopiata, riprodotta, tradotta, distribuita o ridotta a formato leggibile da macchina o mezzo elettronico senza il consenso scritto preventivo di **CTC Analytics AG**, ad eccezione di quanto previsto dalle leggi sul copyright.
PAL is a registered trademark of **CTC Analytics AG** | Switzerland