EVALUATION OF µ-SOLID PHASE EXTRACTION FOR THE ANALYSIS OF PESTICIDE RESIDUES IN CEREALS Elena Hakme^{1;} Mario Francesco Mirabelli²; Ederina Ninga¹; Thomas Preiswerk²; Mette Erecius Poulsen¹ ¹Technical University of Denmark, Copenhagen, Denmark; ²CTC Analytics AG, Zwingen, Switzerland

ABSTRACT

Laboratory automation in routine food control laboratories is increasing with the increased need for cost-effective, time saving, and human error reducing techniques. Novel micro-solid-phase extraction (µ-SPE) extraction devices were recently developed and adopted by some laboratories in a fully automate workflow. The system offers several workflows allowing, among others, automatic preparation of calibration standards, automatic sample dilution, and

automatic clean-up of sample extracts. The system with the implemented workflows was evaluated for its applicability in the analysis of pesticide residues in cereals. The µ-SPE clean-up method was applied to acetonitrile extracts obtained from the QuEChERS sample preparation and was developed for the analysis of pesticide residues in cereals by gas chromatography-Orbitrap mass spectrometry (GC-Orbitrap-MS).

LS3, µ-SPE

1000 µL syringue

LS1, Injection

10 µL syringue

25 µL syringue



GC Parameters		MS Parameters	
Carrier gas Column	He Thermo Scientific Trace	MS Transfer Line temperature	280°C
	length × 0.25 mm i.d. ×	Ion source Temperature	280°C
	0.25 µm	Ionization type	EI
Injection mode Injector temperature	PTV 70°C	Acquisiton mode	Scan Mode
Carrier flow	1.2 mL/min	Scan range	50-550 m/z
Injection volume	1 µL	Mass resolution	60 k (FWHM at
$\begin{array}{c} 350 \\ 300 \\ 250 \\ 200 \\ 150 \\ 100 \\ 50 \\ 0 \\ 2 \\ 4 \\ 6 \\ 8 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 18 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 12 \\ 14 \\ 16 \\ 18 \\ 10 \\ 10 \\ 10 \\ 10 \\ 10 \\ 10 \\ 10$	20 22 24 26 28 30 32 34 36 38 40	Lock masses (column bleed m/z)	355.06993 ; C9 281.05114 ; C7 207.03235 ; C5 133.01356 ; C3 73.04680 ; C3H







Time (mm:ss)	Steps		
0:30	Required tool selected		
	Syringe wash: 2 cycles at wash position 1		
01:00	Load sample onto µ-SPE		
	Perform 3 filling strokes		
03:00	Load sample onto µ-SPE cartridge: 300 µ		
01:00	Syringe wash: 2 cycles at wash position 1		
01:00	Required tool selected		
01:00	Syringe wash: 2 cycles at wash position 1		
	Rinse		
01:00	Add 15 µL internal standard		
	Perform 3 filling strokes		
01:00	Add internal standard: 15 µL		
01:00	Required tool selected		
	Syringe wash: 1 cycle at wash position 2		
01:00	Rinse		
00:30	Move to sample at position 1		
01:00	Perform 3 filling strokes		
	Aspirate 1 µL		
	Inject sample		







pesticides in wheat, 142 pesticides in pesticides in rice. Spike recovery values were 70–120% for all pesticides and the repeatability, calculated as the relative standard deviation, was less than 20%. The limits of quantitation achieved were ensuring compliance with the maximum