

MS²field: High-frequency in situ environmental sampling and HRMS in a transportable container

Michael Stravs¹, Heinz Singer¹, Christian Stamm¹, Christoph Ort²,
Reto Bolliger³, Guenter Boehm³, Thomas Moehring⁴

PAL SYSTEM
Ingenious sample handling

eawag
aquatic research

¹ Department of Environmental Chemistry, Eawag, Dübendorf, Switzerland; ² Department of Urban Water Management, Eawag, Dübendorf, Switzerland; ³ CTC Analytics, Zwingen, Switzerland; ⁴ Thermo Fisher Scientific, Bremen, Germany

OVERVIEW


- Automated micropollutant measurement with filtration – online SPE – LC – HRMS from raw wastewater
- 1 week continuous uptime, 20 min measurement frequency
- Target quantification to ng/L and non-target feature extraction using spectral analysis

INTRODUCTION

Motivation

- Traditional sampling methods** (grab samples, time-proportional samples) require enormous efforts to **capture real-life concentration dynamics** occurring in the environment, e.g.
 - daily/weekly dynamics of pharmaceuticals/PCPs in wastewater influent
 - micropollutant spills in rivers
 - herbicide runoff during rain events
- Manual **sampling and sample preparation** is labor-intensive and error-prone
- Sample **transport and storage** presents major logistical issues

Schematic time profile depending on sampling method



Objectives

- Development of a **transportable automated sampling and measurement platform** to capture micropollutant temporal dynamics
- Automated sampling of **raw wastewater and surface water**
- High-resolution mass spectrometry** to capture a broad compound spectrum and to enable non-target analysis
- High **temporal resolution** (20 min measurement frequency)
- Unsupervised** measurement for **extended periods** (up to 1 week)

Challenges

- Particulate matter: avoid **clogging** etc.
- Robustness**: long unattended operation, error tolerance

Development set-up


- Small-scale experimental **wastewater treatment plant** (200 pers. eq.), influent is sampled after sand trap

Acknowledgements: Richard Fankhauser, Marco Kipf (Eawag); Dr. Nicole Zehethofer, Dr. Siegrun Mohring (Thermo Fisher Scientific)



METHOD

Setup using online SPE – LC – HRMS



Sample preparation:

- Pre-rinse inlet (solenoid ON 30 sec)
- Wash syringe with sample (3 mL)
- Load internal standard to the loop
- Dilute sample (0.75 mL) and load loop

After enrichment:

- Wash loop and sample inlet
- Post-rinse inlet (solenoid ON 90 sec)

Materials:

- SPE: XBridge 10 μ m, 2.0 x 10 mm (Waters)
- LC: XBridge 3.5 μ m, 2.1 x 50 mm (Waters)

Mass spectrometry:


- Polarity switching
- MS1 – top 3/top2 ddMS2
- MS1 – 2x2-window DIA




RESULTS

System evaluation

Time series: Candesartan (antihypertensive drug)




System response time




- LOQs** for 59 micropollutants in wastewater: 19x <10 ng/L, 20x 10-20 ng/L
- Carryover** mostly <1%, for individual compounds up to 10%
- MS sensitivity loss** (spray clogging) after **ca. 1 week** of continuous operation
- Bypass/filter operation max. 1 week** between servicing

Non-target feature extraction




Online live monitoring



- Sensor input** for instrument health
- Automated **data evaluation** after each sample
- MQTT** dashboard allows various data sources and easy visualization

OUTLOOK

In development: SPME – DBDI – HRMS



Deployment in air-conditioned trailer

- Field studies in sewers and rivers

