

A close-up photograph of a laboratory instrument, likely a Solid Phase Microextraction (SPME) system. The image shows various metal components, including a cylindrical part with a 'Tool' label and a vertical component with 'SPME Arrow' and 'Tool' labels. The background is blurred, showing other laboratory equipment.

PFAS in Cosmetics

A Focus on Solid Phase Microextraction (SPME)
Part 2 of the Webinar Series

 **PAL** SYSTEM
Ingenious sample handling

Yes, you are in webinar!



Microphone and webcam are **turned off**




Chat function is **turned off**



Additional information can be found in the Q&A window



Recording in progress



Disclaimer

Please be advised that this webinar is being recorded and later available at the content hub of palsystem.com.

While your microphone and camera will be disabled, your name and any questions you may ask in the Q&A window will be visible to other participants.

By participating in the webinar, you consent to this recording and its use in the content hub.


A webinar series around PFAS – Part 2



Hagen M. Gegner, PhD

Scientific Communications and Projects, CTC Analytics AG

- Background in clinical metabolomics
- Method development in LC-MS and study design
- Focused on sample preparation and quality control



Agenda

Why Cosmetics?

A Direct Exposure Route & An Analytical Challenge

What is SPME?

Principles of a Solvent-Free Microextraction

How is it automated?

A "Walk-Away" Workflow for PFAS in Water

How does it perform?

SPME and PFAS

A Mascara Case Study: SPME vs. μ SPE

Conclusion & Q&A



Ubiquitous PFAS Contamination: From Industry to Your Plate

Di Giorgi, A., et al. (2023). Analysis of perfluoroalkyl substances (PFAS) in conventional and unconventional matrices: Clinical outcomes. Journal of Pharmaceutical and Biomedical Analysis Open, 1, 100002.

Dhiman, S., & Ansari, N. G. (2024). A review on extraction, analytical and rapid detection techniques of Per/Poly fluoro alkyl substances in different matrices. Microchemical Journal, 196, 109667.




PFAS

(Per- and polyfluoroalkyl substances) are a class of over 4,700 synthetic chemicals used since the 1940s.



- Teflon
- Textiles
- Packaging
- Paint

The strength of the Carbon-Fluorine (C-F) bond makes them extremely resistant to heat, chemical, and biological degradation.



A Direct Exposure Route

Di Giorgi, A., et al. (2023). Analysis of perfluoroalkyl substances (PFAS) in conventional and unconventional matrices: Clinical outcomes. Journal of Pharmaceutical and Biomedical Analysis Open, 1, 100002.

Dhiman, S., & Ansari, N. G. (2024). A review on extraction, analytical and rapid detection techniques of Per/Poly fluoro alkyl substances in different matrices. Microchemical Journal, 196, 109667.

Oddný Ragnarsdóttir, Mohamed Abou-Elwafa Abdallah, Stuart Harrad, Dermal bioavailability of perfluoroalkyl substances using in vitro 3D human skin equivalent models, Environment International, Volume 188, 2024, 108772, ISSN 0160-4120, <https://doi.org/10.1016/j.envint.2024.108772>.


A Direct Exposure Pathway

The application of these products represents a direct route for human exposure.



Intentionally Added Ingredients

Unlike many environmental samples where PFAS are unintended contaminants, they are often added to cosmetics on purpose



A Direct Exposure Route

Di Giorgi, A., et al. (2023). Analysis of perfluoroalkyl substances (PFAS) in conventional and unconventional matrices: Clinical outcomes. Journal of Pharmaceutical and Biomedical Analysis Open, 1, 100002.

Dhiman, S., & Ansari, N. G. (2024). A review on extraction, analytical and rapid detection techniques of Per/Poly fluoro alkyl substances in different matrices. Microchemical Journal, 196, 109667.

Oddný Ragnarsdóttir, Mohamed Abou-Elwafa Abdallah, Stuart Harrad, Dermal bioavailability of perfluoroalkyl substances using in vitro 3D human skin equivalent models, Environment International, Volume 188, 2024, 108772, ISSN 0160-4120, <https://doi.org/10.1016/j.envint.2024.108772>.



Cosmetics are a proven exposure route, as a 2024 study confirmed PFAS are absorbed through the skin.

- Shorter-chain PFAS were readily absorbed, with up to 59% of the applied dose permeating the skin barrier.
- Longer-chain PFAS were retained at high concentrations within the skin tissue, even if not directly absorbed.

A Direct Exposure Pathway
The application of these products represents a direct route for human exposure.



Intentionally Added Ingredients

Unlike many environmental samples where PFAS are unintended contaminants, they are often added to cosmetics on purpose



The Analytical Challenge of Mascara

Why is Analysis So Difficult?

Di Giorgi, A., et al. (2023). Analysis of perfluoroalkyl substances (PFAS) in conventional and unconventional matrices: Clinical outcomes. *Journal of Pharmaceutical and Biomedical Analysis Open*, 1, 100002.

Dhiman, S., & Ansari, N. G. (2024). A review on extraction, analytical and rapid detection techniques of Per/Poly fluoro alkyl substances in different matrices. *Microchemical Journal*, 196, 109667.

Oddný Ragnarsdóttir, Mohamed Abou-Elwafa Abdallah, Stuart Harrad, Dermal bioavailability of perfluoroalkyl substances using in vitro 3D human skin equivalent models, *Environment International*, Volume 188, 2024, 108772, ISSN 0160-4120, <https://doi.org/10.1016/j.envint.2024.108772>.




Trace Concentrations

PFAS are often present at very low levels (ng/g or parts-per-billion), demanding highly sensitive analytical instruments.



Matrix Effects

During analysis, co-extracted ingredients can artificially suppress or enhance the PFAS signal in the mass spectrometer, a common challenge requiring careful method design.



Matrix Complexity

Cosmetics are a dense mixture of waxes, oils, pigments, and polymers that can interfere with extraction and lead to inaccurate results.



Background Contamination

The ubiquity of PFAS in laboratory environments (e.g., from PTFE tubing) requires strict contamination controls and specialized equipment to ensure reliable results.



The Solution: A Deeper Look at SPME

What is Solid Phase Microextraction (SPME)?

Jiang, G., Huang, M., Cai, Y., Lv, J., & Zhao, Z. (2006). Progress of Solid-Phase Microextraction Coatings and Coating Techniques. Journal of Chromatographic Science, 44(6), 324–333.
Dugheri, S., Mucci, N., Cappelli, G., Trevisani, L., Bonari, A., Bucaletti, E., Squillaci, D., & Arcangeli, G. (2022). Advanced Solid-Phase Microextraction Techniques and Related Automation: A Review of Commercially Available Technologies. Journal of Analytical Methods in Chemistry, 2022, 8690569.

The Principle


SPME is a solvent-free sample preparation technique developed by Prof. Janusz Pawliszyn (1989). It works on the principle of equilibrium partitioning.



SPME ARROW

Extraction

The way analytes are extracted depends on the fiber's coating material, which can be either a liquid-like polymer or a solid porous layer.



Adsorption

This is the primary mechanism for solid, porous coatings like Divinylbenzene (DVB) or Carboxen (CAR). Analytes accumulate on the outer surface of the sorbent material



Absorption

Analytes partition from the sample and diffuse into the entire volume of the coating. This occurs with liquid-like coatings such as polydimethylsiloxane (PDMS) or polyacrylate (PA).


The Solution: A Deeper Look at SPME

What is Solid Phase Microextraction (SPME)?

Jiang, G., Huang, M., Cai, Y., Lv, J., & Zhao, Z. (2006). Progress of Solid-Phase Microextraction Coatings and Coating Techniques. *Journal of Chromatographic Science*, 44(6), 324–333.
 Dugheri, S., Mucci, N., Cappelli, G., Trevisani, L., Bonari, A., Bucaletti, E., Squillaci, D., & Arcangeli, G. (2022). Advanced Solid-Phase Microextraction Techniques and Related Automation: A Review of Commercially Available Technologies. *Journal of Analytical Methods in Chemistry*, 2022, 8690569.


The Principle

SPME is a solvent-free sample preparation technique developed by Prof. Janusz Pawliszyn (1989). It works on the principle of equilibrium partitioning.




Key Differences

- The Arrow is fabricated with a robust stainless-steel core and a protective outer sheath. This design confers superior mechanical stability and durability, reducing the risk of fracture common to the more fragile fused-silica core of standard fibers.
- The Arrow possesses a significantly greater sorptive phase volume than traditional fibers. This increased volume provides a higher analyte loading capacity, which can lead to a greater preconcentration factor from the sample matrix.




The Automated Workflow in 3 Key Steps



Conditioning & Sample Preparation


The SPME fiber is first heated in a conditioning station to ensure it is clean and ready for extraction, preventing sample-to-sample carryover.

Simultaneously, the autosampler can prepare the next sample, for example by moving it to the agitator.



The Automated SPME Workflow Principle

The Automated Workflow in 3 Key Steps



Conditioning & Sample Preparation



The SPME fiber is first heated in a conditioning station to ensure it is clean and ready for extraction, preventing sample-to-sample carryover.

Simultaneously, the autosampler can prepare the next sample, for example by moving it to the agitator.


Incubation & Extraction

The sample vial is heated and agitated for a set time to help the PFAS partition out of the complex matrix.

The conditioned SPME is then automatically exposed to the sample (either via direct immersion or headspace) for a precisely controlled extraction time.




The Automated Workflow in 3 Key Steps



Conditioning & Sample Preparation

The SPME fiber is first heated in a conditioning station to ensure it is clean and ready for extraction, preventing sample-to-sample carryover.


Simultaneously, the autosampler can prepare the next sample, for example by moving it to the agitator.



Incubation & Extraction

The sample vial is heated and agitated for a set time to help the PFAS partition out of the complex matrix.

The conditioned SPME is then automatically exposed to the sample (either via direct immersion or headspace) for a precisely controlled extraction time.




Desorption & Injection

The SPME, now loaded with the concentrated PFAS, is withdrawn from the sample and ready for analysis.


For GC-MS: Desorption is a single, rapid thermal step directly inside the hot GC injector.

For LC-MS: Desorption is a two-step process: first, the analytes are desorbed using a solvent in a separate vial, and then that solvent is injected into the instrument



The Automated SPME Workflow Application

Optimizing the Extraction: Time is Key




Extraction Time Profile for Volatile PFAS using Direct Immersion SPME (DI-SPME).

This graph shows the amount of different volatile PFAS extracted over 90 minutes, measured as GC-MS peak area. For this experiment, a DVB/Car/PDMS SPME fiber was immersed directly into ultrapure water spiked with the target analytes. The extraction was performed at 60 °C with an agitation of 250 rpm. This type of profile is critical for method development, as it helps determine the optimal time needed to extract enough analyte for sensitive detection. Under these conditions, most of the tested PFAS approach equilibrium after approximately 60 minutes.



Automation is key

Its precise control of all parameters allows for highly reproducible results and efficient overlapping of sample prep and injection.





Extraction is a dynamic process

The amount of analyte captured by the SPME increases over time until it reaches equilibrium with the sample.

- For Direct Immersion (DI-SPME), equilibrium for many volatile PFAS is reached in approximately 60 minutes.
- For Headspace (HS-SPME), the process is significantly faster, with most analytes reaching equilibrium in about 25 minutes.



Optimizing the Extraction: The Effect of Temperature



The Effect of Extraction Temperature on SPME Efficiency for Volatile PFAS.

These charts compare the extraction efficiency (measured as GC-MS peak area) for volatile PFAS at temperatures from 30 °C to 60 °C using two different modes: **(A) Direct Immersion (DI-SPME)** and **(B) Headspace (HS-SPME)**.


In Direct Immersion (A), extraction efficiency increases with temperature for all compounds, which is particularly evident for the less volatile EtFOSA and MeFOSE.

In Headspace (B), the trend is different; while extraction of EtFOSA and MeFOSE still increases with temperature, the extraction of the more volatile FTOHs peaks at 50 °C and then declines. This illustrates the need to carefully optimize temperature based on the specific analytes and the chosen extraction mode.

Account for contamination sources

There are upgrade kits available for PFAS Analysis. They focus on the tubing, LC/MS Tool, vials and syringes.

Always run blanks for your QC/QA strategies.



Temperature directly influences the efficiency of the SPME process and must be carefully optimized for target analytes.

- In Direct Immersion (DI-SPME), higher temperatures generally increase the amount of analyte extracted by speeding up diffusion, especially for less volatile compounds. An optimal temperature (e.g., 60 °C) is often chosen to maximize recovery while maintaining good reproducibility.
- In Headspace (HS-SPME), the effect is more complex, creating a trade-off: Higher temperatures help less volatile compounds move into the headspace, increasing their extraction. However, because sorption onto the fiber is an exothermic process, temperatures that are too high can reduce the extraction of the most volatile compounds.

The Automated SPME Workflow Application

Automated Analysis of Volatiles in Water





Table 4: Summary of PFAS calibration range and linearity results

Compound	Calibration Range (ng/L)	R ²	RF (Response Factor) %RSD
PFHxI	2.5-2000	0.993	10.89
PFOI	2.5-2000	0.997	10.26
4:2 FTI	2.5-800	0.993	8.28
6:2 FTI	25-800	0.994	13.53
8:2 FTOH	25-2000	0.997	5.37
6:2 FTAC	25-2000	0.998	19.87
8:2 FTI	2.5-800	0.996	13.59
10:2 FTOH	2.5-2000	0.999	10.38
6:2 FTMAC	2.5-800	0.995	12.43
8:2 FTAC	5-250	0.995	14.81
8:2 FTMAC	2.5-250	0.998	19.51
MeFOSA	5-2000	>0.999	17.79
EtFOSA	10-2000	0.999	11.40

A fully automated and cost-effective method

The selected extraction temperature (50 °C) and time (30 min) represent an effective compromise, providing good sensitivity across a diverse group of PFAS with varying volatilities in a single, unattended run.



Sample Preparation

Water samples were "salted out" with 2% NaCl to enhance the extraction efficiency of the target compounds.

Automated Incubation & Extraction

The autosampler incubated the sample for 5 minutes at 50 °C and then exposed the SPME fiber to the vial's headspace for 30 minutes.

Desorption & Analysis

The fiber was transferred directly to the hot GC inlet (240 °C) for a rapid 7-minute thermal desorption.


Analysis of Volatile PFAS in Water Using Head-Space Solid Phase Microextraction-Gas Chromatography/Mass Spectrometry (HS-SPME GC/MS) Source: Shimadzu Scientific Instruments



Tackling the Complex Mascara Matrix

Choosing the Right Tool for the Job

SPME Fiber Workflow



SPME Fiber Workflow


This workflow was performed manually in the study to establish a proof-of-concept for this challenging matrix.

- 0.5 g of mascara was dispersed in 5 mL of 100% water.
- The SPME fiber was immersed directly into the mascara-water mixture for 60 minutes to extract the PFAS.
- Analytes were desorbed from the fiber using a small amount of solvent, which was then collected for analysis.



Olomukoro, A. A., Lüthy, L., Flug, T., & Gionfriddo, E. (2025). Evaluation of extraction methodologies for PFAS analysis in mascara: a comparative study of SPME and automated μ SPE. Analytical and Bioanalytical Chemistry. <https://doi.org/10.1007/s00216-025-05908-x>

μ SPE Workflow



μ SPE Workflow

The core extraction and elution steps of this workflow were automated on a PAL System

- 0.5 g of mascara was dispersed in 5 mL of 100% methanol to dissolve the sample matrix.
- The mascara-methanol mixture required manual filtration before being placed on the autosampler to prevent the μ SPE cartridges from clogging.
- The PAL System then handled the conditioning of the μ SPE cartridge, loading of the filtered sample, and elution of the final extract.

Both workflows benefit significantly from automation.

SPME Fibers/Arrows can be fully automated for unattended batch processing. The core steps of μ SPE are automated, and the initial filtration step can also be integrated into the workflow for a complete walk-away solution.



Choosing the Right Tool for the Job



SPME (Fiber)

Extraction

Equilibrium Partitioning

Sample Pre-treatment

No filtration required

Optimal Dispersion Solvent

100% Water

Sensitivity (LOQ)

Better for hydrophilic PFAS (e.g., PFBS)

Device Reusability

Reusable Fiber



Micro-SPE Cartridges

Extraction

Exhaustive Extraction (SPE)

Sample Pre-treatment

Filtration of sample is required

Optimal Dispersion Solvent

100% Methanol

Sensitivity (LOQ)

Better for hydrophobic PFAS (e.g., PFOS)

Device Reusability

Single-use Cartridges

Both methods are effective for real-world product safety testing, successfully detecting and quantifying regulated PFAS in commercial mascara products.

- SPME offers a simpler, filtration-free workflow, while the automated μ SPE method required a manual filtration step to handle the complex mascara matrix. The techniques are complementary, not competing.
- SPME showed superior sensitivity for hydrophilic PFAS, whereas automated μ SPE was more effective for key hydrophobic PFAS.

"Analytical Toolbox" concept, where the best strategy is to choose the right combination of tools for the job. For this matrix, SPME had a simpler workflow for hydrophilic targets, while automated μ SPE offered higher sensitivity for hydrophobic ones. For instrumentation, LC-MS is the standard for ionizable PFAS, while GC-MS excels for the neutral and volatile species.



Olomukoro, A. A., Lüthy, L., Flug, T., & Gionfriddo, E. (2025). Evaluation of extraction methodologies for PFAS analysis in mascara: a comparative study of SPME and automated μ SPE. Analytical and Bioanalytical Chemistry. <https://doi.org/10.1007/s00216-025-05908-x>

Both workflows benefit significantly from automation.

SPME Fibers/Arrows can be fully automated for unattended batch processing. The core steps of μ SPE are automated, and the initial filtration step can also be integrated into the workflow for a complete walk-away solution.



SPME is a Powerful and Versatile Tool



Cosmetics are a Direct Exposure Route

PFAS are absorbed through the skin. Studies show shorter-chain PFAS can be readily absorbed, while longer-chain compounds are retained within the skin tissue itself.




SPME is a Powerful and Versatile Tool

Solid Phase Microextraction is a sensitive, solvent-free technique that effectively extracts PFAS from complex matrices, making it an excellent choice for modern, environmentally conscious labs.



SPME performance depends on the **precise optimization** of factors such as extraction time and temperature. Automation secures these ideal conditions for each sample, guaranteeing exceptional reproducibility and facilitating dependable, high-throughput "walk-away" operation.



SPME is a Powerful and Versatile Tool



Cosmetics are a Direct Exposure Route

PFAS are absorbed through the skin. Studies show shorter-chain PFAS can be readily absorbed, while longer-chain compounds are retained within the skin tissue itself.



SPME is a Powerful and Versatile Tool

Solid Phase Microextraction is a sensitive, solvent-free technique that effectively extracts PFAS from complex matrices, making it an excellent choice for modern, environmentally conscious labs.



SPME performance depends on the **precise optimization** of factors such as extraction time and temperature. Automation secures these ideal conditions for each sample, guaranteeing exceptional reproducibility and facilitating dependable, high-throughput "walk-away" operation.





Build the Right Analytical Toolbox

There is no single best method for all PFAS; the most effective strategy is to choose the right combination of tools for the job.

- For instrumentation, LC-MS is the standard for ionizable PFAS, while GC-MS excels for the neutral and volatile species that LC can miss.
- For sample preparation, the mascara study showed that SPME offers a simpler, filtration-free workflow ideal for hydrophilic PFAS, while automated μ SPE provides higher sensitivity for hydrophobic ones.

Importantly, all sample preparation aspects can be automated.





Ask your question in the **Q & A window** right now



Use the form linked in the pinned message



Send us an email **info@palsystem.com**



Scientific References

1. Di Giorgi, A., et al. (2023). *Journal of Pharmaceutical and Biomedical Analysis Open*, 1, 100002.
2. Dhiman, S., & Ansari, N. G. (2024). *Microchemical Journal*, 196, 109667.
3. Dugheri, S., et al. (2022). *Journal of Analytical Methods in Chemistry*, 2022, 8690569.
4. Jiang, G., et al. (2006). *Journal of Chromatographic Science*, 44(6), 324-333.
5. Lehotay, S. J., et al. (2016). *Chromatographia*, 79, 1113-30.
6. Martínez-Pérez-Cejusia, H., et al. (2025). *Analytica Chimica Acta*, 1345, 343746.
7. Olomukoro, A. A., et al. (2025). *Analytical and Bioanalytical Chemistry*, <https://doi.org/10.1007/s00216-025-05908-x>
8. Ragnarsdóttir, O., et al. (2024). *Environment International*, 188, 108772.
9. Shimadzu Scientific Instruments, Application Note GCMS-2402, "Analysis of Volatile PFAS in Water Using HS-SPME-GC/MS".

PAL System Resources

- PAL SPME Arrow: [Discover the Power of SPME and the Benefits of SPME Arrows](#)
- PAL μ SPE: [Micro-SPE - A Green Extraction technique in modern laboratories](#)
- Upgrade Kits for PFAS Analysis: [PFAS Analysis - Equip your PAL System for PFAS Analysis](#)
- PAL System Content Hub: [PAL Content Hub - Search our Library](#)