



PAL System Vials & Vial Caps



- Borosilicate hydrolytic glass vials and polypropylene vials available
- Better sealing due to flat finish of contact surface to septum
- Improved bottom shape and thicker glass walls for fast and reliable transportation
- Reduced glass breakage through refined high quality glass
- Large labelling area for chemical resistant pencil inscriptions
- PAL System tested and approved

Please see below all our available Vials and Vial Caps:

Vials	
PN CTC	Description
Vial-0.3-ND9-PP-100	PAL System Micro-Vial 2CV, 0.3mL Clear PP, designed for the PAL Autosampler. 12x32mm, with conical insert, fits ND9 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-1.5-ND9-AG-100*	PAL System Vial 2CV, 1.5ml Amber Glass with Label, designed for the PAL Autosampler. 12x32mm, 1st Class Hydrolytic Glass, fits ND9 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-1.5-ND9-CG-100*	PAL System Vial 2CV, 1.5ml Clear Glass with Label, designed for the PAL Autosampler. 12x32mm, 1st Class Hydrolytic Glass, fits ND9 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-1.5-ND9-PP-100	PAL System Vial 2CV, 1.5mL Clear PP, designed for the PAL Autosampler. 12x32mm, with filling lines, fits ND9 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-10-DC20-CG-100	PAL System Vial 10CV, 10ml Clear Glass with Label, for PAL Wash Modules. 46x22.5mm, 1st Class Hydrolytic Glass, DIN-Crimp Neck, Pk of 100 Pcs, clip caps not included
Vial-10-ND18-AG-100*	PAL System Vial 10CV, 10ml Amber Glass with Label, designed for the PAL Autosampler. 46x22.5mm, 1st Class Hydrolytic Glass, with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-10-ND18-CG-100*	PAL System Vial 10CV, 10ml Clear Glass with Label, designed for the PAL Autosampler. 46x22.5mm, 1st Class Hydrolytic Glass, with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-20-DC20-AG-100	PAL System Vial 20CV, 20ml Clear Glass for Headspace. 75.5x22.5mm, 1st Class Hydrolytic Glass, DIN-Crimp Neck, Pk of 100pcs, screw caps not included
Vial-20-DC20-CG-100	PAL System Vial 20CV, 20ml Clear Glass for Headspace. 75.5x22.5mm, 1st Class Hydrolytic Glass, DIN-Crimp Neck, Pk of 100pcs, screw caps not included
Vial-20-ND18-AG-100*	PAL System Vial 20CV, 20ml Amber Glass with Label, designed for the PAL Autosampler. 75.5x22.5mm, 1st Class Hydrolytic Glass with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-20-ND18-CG-100*	PAL System Vial 20CV, 20ml Clear Glass with Label, designed for the PAL Autosampler. 75.5x22.5mm, 1st Class Hydrolytic Glass with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs, screw caps not included
Vial-40-ND24-CG-100	PAL System Vial 40CV, 40ml Clear Glass with Label, designed for the PAL Autosampler. 95x27.5mm, 1st Class Hydrolytic Glass with flat top for optimized sealing. For ND24 Screw Caps, Pk of 100 Pcs, screw caps not included
Caps	
PN CTC	Description
Cap-ND9-St-SP10Sb-100	PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Magnetic, Gold, Silicone/PTFE/Starburst-Slitted Septa 1.0mm, Pk of 100 pcs
Cap-ND9-St-SP10-100*	PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Magnetic, Gold, Silicone/PTFE Septa 1.0mm Pk of 100 pcs
Cap-ND9-PP-SP10S-100	PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Silicone/PTFE/I-Slitted Septa 1.0mm, Pk of 100 pcs
Cap-ND9-PP-SP10-100	PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Silicone/PTFE Septa 1.0mm, Pk of 100 pcs
Cap-DC20-PE8-SP13Y-100	PAL System Clip Cap for PAL Wash Modules. 22mm PE with 8mm center hole, Silicone/PTFE/Y-Slitted Septa 1.3mm Pk of 100 pcs
Cap-ND18-StC-SP20-100	PAL System Screw Re-Cap 10CV and 20CV, designed for the PAL Autosampler. ND18 Magnetic, PTFE Metal Coating for optimized screw and sealing performance, Silicone/PTFE Septa 2.0mm, Pk of 100 pcs
Cap-ND18-St-SP15-100*	PAL System Screw Cap 10CV and 20CV, designed for the PAL Autosampler. ND18 Magnetic, Silver, Silicone white/PTFE blue Septa 1.5mm, Pk of 100 pcs
Cap-ND24-PP-SP32-100	PAL System Screw Cap 40CV, designed for the PAL Autosampler. ND24, Silicone/PTFE Septa 3.2mm, Pk of 100 pcs
Trial Set:	
Part No.	Description
PAL3-Vial-Sample-Set	Vial and Cap Sample Set consisting of: 2pc 1.5ml transparent clear glass and 2pc 1.5ml amber glass vial, fits ND9 Screw Cap 2pc 10ml transparent clear glass and 2pc 10ml amber glass vial, fits ND18 Screw Cap 2pc 20ml transparent clear glass and 2pc 20ml amber glass vial, fits ND18 Screw Cap 4pc Screw Cap, ND9, Magnetic, Slitted Septum 1.0mm 4pc Screw Cap, ND18, Magnetic 4pc Screw Re-Cap, ND18, Magnetic

© 2020 CTC Analytics AG. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine-readable form without the prior written permission from CTC Analytics AG, except as permitted under copyright laws.

CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

PAL is a registered trademark of CTC Analytics AG | Switzerland

Date of print: 04.2020 | CTC Analytics AG, Industriestrasse 20, CH-4222 Zwingen, Switzerland, www.palsystem.com, Contact: info@ctc.ch