

# A single-sample workflow for joint metabolomic and proteomic analysis of clinical specimens

Hagen M. Gegner<sup>1,7</sup>, Thomas Naake<sup>2</sup>, Karim Aljakouch<sup>3</sup>, Aurelien Dugourd<sup>4</sup>, Georg Kliewer<sup>3</sup>, Torsten Müller<sup>3</sup>, Dustin Schilling<sup>4</sup>, Marc A. Schneider<sup>5</sup>, Nina Kunze Rohrbach<sup>1</sup>, Thomas G.P. Grünewald<sup>6</sup>, Rüdiger Hell<sup>1</sup>, Julio Saez Rodriguez<sup>4</sup>, Wolfgang Huber<sup>2</sup>, Gernot Poschet<sup>1</sup> & Jeroen Krijgsveld<sup>3</sup>



## Overview

Integrating proteomic and metabolomic data enhances our understanding of biological processes and disease mechanisms. Traditional multi-omics workflows often require separate sample processing, leading to increased variability. To address this, we developed MTBE-SP3, an optimized single-sample workflow that combines a 75% EtOH/MTBE metabolite extraction with automated single-pot solid-phase-enhanced sample preparation for proteins. **This approach reduces variability and improves robustness in multi-omics analyses**, allowing for deeper insights into cellular regulation and disease pathways. Here, we present the MTBE-SP3 workflow and its effectiveness using lung adenocarcinoma patient samples, along with the potential for further automation to boost throughput and standardization in multi-omics research.

## Key benefits of MTBE-SP3



- Single-sample processing for **reduced variability** and **more coherent** multi-omics analysis.
- Broadly tested applicability to various biological matrices like FFPE, fresh-frozen tissues, plasma, serum and cells as well as a pilot study.
- Potential for full automation that minimizes technical variability for clinical workflows.

## Pilot Study

### Lung cancer



Applied to lung adenocarcinoma patients (10) we identified distinct protein and metabolite changes between tumor (TT) and normal tissues (NT). **Integrating these molecular layers revealed additional dependencies** and emphasized mitochondrial dysfunction in tumor development, demonstrating the value of multi-omics analyses. Finally, we showed that single sample integration results in **better correlation** in the ocEAn analysis.



## Single Sample Network Analysis

- Integration revealed **additional pathways** that are affected in Tumor Tissue vs Normal Tissue
- Network Analysis shows better correlation for single sample data

autoSP3  
Protein data showed clear differences between Tumor Tissue vs Normal Tissue



## Why EtOH/MTBE?

**Lower variance** than Bligh & Dyer  
Broad metabolite coverage  
Protein pellet at the bottom for **automation**

## What is the overlap in percentage? AutoSP3 vs MTBE-SP3


- 85%** FFPE
- 89%** Fresh-frozen
- 91%** Plasma
- 98%** Cells



See the extraction comparison online.

## Next Step: More Automation

Automated adjusted **Matyash extraction (MTBE) with 75% EtOH**



Explore more automation options



See the open access publication

**Affiliations**  
 1. Centre for Organismal Studies (COS), Metabolomics Core Technology Platform, Heidelberg University, Im Neuenheimer Feld 360, 69120, Heidelberg, Germany  
 2. Genome Biology Unit, European Molecular Biology Laboratory (EMBL), Meyerhofstr. 1, 69117, Heidelberg, Germany  
 3. Faculty of Medicine, Heidelberg University, 69120, Heidelberg, Germany and Division Proteomics of Stem Cells and Cancer, German Cancer Research Center (DKFZ), Im Neuenheimer Feld 581, 69120, Heidelberg, Germany  
 4. Institute for Computational Biomedicine, Bioquant, Faculty of Medicine, Heidelberg University and Heidelberg University Hospital, Im Neuenheimer Feld 130, 69120, Heidelberg, Germany  
 5. Translational Research Unit, Thoraxklinik at Heidelberg University Hospital, Röntgenstraße 1, 69126, Heidelberg, Germany and Translational Research Center Heidelberg (TRC), Member of The German Center for Lung Research (DZL), Im Neuenheimer Feld 156, 69120, Heidelberg, Germany  
 6. Division of Translational Pediatric Sarcoma Research, German Cancer Research Center (DKFZ), German Cancer Consortium (DKTK), Im Neuenheimer Feld 280, 69120, Heidelberg, Germany and Hopp-Children's Cancer Center (KITZ), Heidelberg, Germany and Institute of Pathology, Heidelberg University Hospital, Im Neuenheimer Feld 224, 69120, Heidelberg, Germany and National Center for Tumor Diseases (NCT), NCT Heidelberg, DKFZ and Heidelberg University Hospital, Heidelberg, Germany  
 7. CTC Analytics AG, Zwingen, Switzerland